

## TENStem eco basic

### PAIN MANAGEMENT IN HOME TREATMENT

The TENStem eco basic is the most popular TENS unit in the schwa-medico range for price sensitive markets. The great diversity of programs offers an easy selection of the treatment parameters. The TENStem eco basic includes 12 pre-set programs such as burst, frequency modulation, muscle stimulation, dynamic stimulation, HAN stimulation.

### 4 FREE PROGRAMS

The TENStem eco basic also offers 4 free programs for own settings. The latest program P16 now offers the possibility of stimulating with monophasic current. P16 is suitable for iontophoresis treatment. TENStem eco basic is ideal for home treatment pain management. It is a valuable TENS unit which offers many features at a very reasonable price.

### SUPPLIED WITH THE FOLLOWING ACCESSORIES

104741	Cable type 7
283400	Self-adhesive electrodes 50 x 50 mm
604000	1.5 V micro battery
	Carton box
100863	Operating manual

### SEE OPTIONAL ACCESSORIES ON PAGES 26 TO 29



### FEATURES



### TYPE OF PROGRAM

<input checked="" type="checkbox"/> PAIN	<input checked="" type="checkbox"/> MUSCLE STIM	<input checked="" type="checkbox"/> DYNAMIC
<input checked="" type="checkbox"/> HAN	<input type="checkbox"/> DENTAL	<input type="checkbox"/> UROLOGICAL
<input type="checkbox"/> SPORT	<input type="checkbox"/> BEAUTY	<input type="checkbox"/> WELLNESS
<input type="checkbox"/> REHABILITATION	<input checked="" type="checkbox"/> RELAXATION	<input type="checkbox"/> SPECIAL

### TECHNICAL DATA

Item no.	104041
Output current	60 mA (with 1 kΩ real)
Frequency	1 – 120 Hz
Pulse width	75 – 300 μs
Power supply	4x AAA batteries